ES INSTRUCCIONES DE USO	INSTRUCCIONES DE USO DE IMPLANTES
EN INSTRUCTIONS FOR USE	IMPLANT INSTRUCTIONS FOR USE
DE GEBRAUCHSANWEISUNG	GEBRAUCHSANWEISUNG ZU DEN IMPLANTATEN
FR INSTRUCTIONS D'UTILISATION	INSTRUCTIONS D'USAGE DES IMPLANTS
IT ISTRUZIONI PER L'USO DEGLI IMPIANTI	ISTRUZIONI PER L'USO DEGLI IMPIANTI
PT INSTRUÇÕES DE UTILIZAÇÃO	INSTRUÇÕES DE UTILIZAÇÃO DE IMPLANTES
TR KULLANIM KILAVUZU	İMLANT KULLANIM KILAVUZU

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Instructions for use applicable to Ziacom® dental implant systems

IMPLANT SYSTEMS		
BRAND	MODEL	MODEL DESCRIPTION
ZIACOM®	ZICI	ZINIC
	ZECI	ZM4
	Z1EC	ZM1
	Z4EC	ZM4 MT
	ZMIC	ZINIC MT
	ZSIC	ZINIC SX
	ZXIC	ZINIC MTX
	Z8EC	ZM8
	ZGCC	GALAXY
	ZVCC	ZV2
	Z1CC	ZM10
	ZKCI	ZMK
	ZMCI	ZMR

IMPORTANT INFORMATION

Carefully read these instructions for use before using the product.

This document contains basic information for the use of original Ziacom® dental implant systems, hereinafter referred to as “Ziacom® dental implants” or simply “Ziacom® products”. This documentation has been prepared as a quick reference guide for the healthcare professional responsible for the treatment, hereinafter referred to as the “User”, and does not provide sufficient indications and technical specifications for the correct use of Ziacom® products. It is neither an alternative nor a substitute for specialized training and professional clinical experience.

Ziacom® products must be used following appropriate treatment planning and strict adherence to the surgical and/or prosthetic protocols established by Ziacom Implants SLU. Carefully read the specific surgical and prosthetic protocols for each product, as well as the instructions for use and maintenance before using any Ziacom® product. These documents may be consulted on our website www.ziacom.com or requested from your official distributor authorized by Ziacom Implants SLU.

1. System Description

Ziacom® products consist of dental implants, prosthetic abutments or prosthodontic components, and surgical and prosthetic instruments.

Ziacom Implants SLU has developed various dental implant systems for performing dental treatments in accordance with the existing scientific literature and current clinical standards. Ziacom Implants SLU develops, manufactures and markets state-of-the-art implant products and a wide range of restorative solutions to meet the specific needs of each patient.

2. Product Description

Ziacom® dental implants are manufactured from biocompatible Grade 4 titanium in accordance with ASTM F67 and ISO 5832-2 standards, except for Zinic MTX and Zinic SX, which are manufactured from Grade 5 ELI titanium in accordance with ASTM F136 and ISO 5832-3 standards. All implants feature a Titansure or Titansure Active surface treatment. Ziacom® dental implants are manufactured in different diameters, lengths, shapes and connections depending on their intended purpose.

3. Indications for Use

Before performing treatment with Ziacom® products, it is essential to evaluate the patient through clinical diagnosis, imaging assessment and model analysis, as well as to determine any possible contraindications to treatment. Successful dental implant treatment must be based on a diagnosis, a treatment plan and a prosthetic rehabilitation design that meets the patient’s needs and expectations.

Ziacom® dental implants are devices intended for the total or partial replacement of teeth in order to restore masticatory function, as well as to provide improved occlusion and distribution of occlusal loads, phonation, prosthetic stability and preservation of the alveolar process. They are supplied in various designs and in different lengths and diameters to accommodate different clinical situations.

Ziacom® dental implants are single-use devices.

Carefully review the clinical indications for each implant model in the corresponding catalogues and surgical protocols before planning surgery. Do not use Ziacom® dental implants in tooth positions or clinical situations that are not described as indications in the corresponding catalogues or surgical protocols.

4. Contraindications

A preoperative medical examination of the patient is required to determine whether absolute and/or relative contraindications or risk factors for dental implant treatment exist. The User is responsible for evaluating the potential benefits and risks of implant treatment and for making the appropriate decisions according to each case. Ziacom® products must not be used in patients who do not meet the minimum medical conditions required for dental implant treatment and prosthetic rehabilitation.

Relative Contraindications

Patient age – Uncontrolled systemic disease – Poor oral hygiene – Active local infection – Pregnancy – Bisphosphonate therapy – Sjögren’s syndrome – Uncontrolled endocrine disease – Head and neck radiation – Stainless steel allergy – Stress – Psychological factors – Bruxism/Parafunction – Prosthetic heart valves – Uncontrolled periodontal disease – Smoking – Alcoholism – Drug addiction – Addictions in general.

Absolute Contraindications

The manufacturer DOES NOT recommend implant treatment in the following cases: Cardiovascular diseases or disorders – Hemophilia – Liver cirrhosis – Severe osteoporosis – Terminal illnesses – Epilepsy – Neoplasms treated with chemotherapy – Blood dyscrasias – Cancer – Uncontrolled diabetes – Untreated periodontitis – Infectious diseases – Leukemia – Immune system disorders – Immunosuppressive treatment – Titanium allergy – Active radiotherapy treatments – Long-term corticosteroid treatments – Anticoagulant treatments – Involuntary muscular disorders – Severe liver disorders – Endocrine disorders – Psychiatric and/or psychological disorders.

Ziacom Implants SLU will update these contraindications when necessary based on information obtained through clinical studies and post-market surveillance.

5. Intended Patients

These medical devices are intended for patients who, in the opinion of the clinician evaluating the case, require treatment in the mandible and/or maxilla. Their intended purpose is the restoration of masticatory function, while also improving occlusion and the distribution of occlusal loads, phonation, aesthetics, prosthetic stability and preservation of the alveolar process or alveolar ridge, together with a positive psychological effect.

Patients must have completed their growth and development. A minimum age of eighteen years is recommended as a guideline; however, cessation of growth may be confirmed by means of a wrist radiograph. Young patients shall be expressly informed of any changes that may need to be made to the prostheses throughout their lifetime, together with the associated aesthetic and economic implications, due to the osseointegration of the implant within living tissue.

6. Storage

Ziacom® products must be stored in their original packaging at a temperature between -10°C and 55°C, in a dry, clean location protected from sunlight and adverse environmental conditions.

7. Packaging and Supply Condition

Ziacom® dental implants are supplied in a cardboard box with a product identification label that serves as a tamper-evident seal. Depending on the reference, the packaging may consist of a blister pack or a container with the following characteristics:

- Heat-sealed plastic blister with removable labels on the back indicating the manufacturer, the complete product reference and the lot number. Inside the blister:

- No Mount models: Implant vial containing the implant
- Mount models: Implant vial containing the implant together with its Mount, clinical screw and cover screw.
- Models with plastic carrier: Implant vial containing the implant attached to the plastic Mount.

- Plastic container containing the implant vial in which the implant is housed and which contains a saline solution for the optimal preservation of its properties.

- The removable labels are intended for implant traceability control, with one label being affixed to the patient’s implant card and the remaining label retained for the clinic’s records.

Ziacom Implants SLU guarantees that all its products undergo a comprehensive manufacturing, inspection and cleaning process before packaging. Before using Ziacom® products, the integrity of the packaging must be verified and it must be ensured that it is not damaged. If any defect or damage is detected, the product must not be used and the issue must be reported to Ziacom Implants SLU or its authorized official distributor.

8. Cleaning and Sterilization

Ziacom® dental implants are supplied sterile, having been sterilized by Beta irradiation at 25 kGy. If the sterile packaging has been altered, tampered with, damaged or opened, sterility can no longer be guaranteed and the product must not be used. Do not use the product after the expiry date indicated on the packaging.

Re-sterilization of Ziacom® dental implants is NOT permitted.

Important

The blister pack or plastic container containing the implant vial is NOT sterile on its outer surface; only the interior of the blister is sterile. Do not place the blister pack within the sterile surgical field. It must NOT be opened until immediately before implant placement.

9. Precautions**Surgery and Instrumentation**

The surgical protocols describe in detail the precautions and important factors to be considered during the surgical placement of dental implants. Prior to the procedure, it is essential to consult the specific surgical protocol for the product, as well as all instructions regarding its correct use contained in the Ziacom® product catalogues.

Particular attention must be paid and appropriate measures taken to prevent instruments, due to their small size, from being swallowed or aspirated by the patient. All manually operated Ziacom® instruments incorporate retention features for connection to each other or for use with dental floss or suture thread in order to prevent accidental ingestion.

It is important to minimise surgical time and damage to the implant recipient site, paying particular attention to bone trauma and maintenance of aseptic conditions within the surgical field. Proper preparation of the surgical site is a fundamental factor for treatment success. This process requires strict adherence to the drilling protocol specified by the manufacturer for each implant type and the correct use of surgical drills with regard to both the number of uses and operating speeds. Surgical drills and taps must be used under abundant and continuous irrigation to prevent bone overheating.

Special care should be taken not to apply excessive force or torque during implant insertion. It is important to monitor implant insertion torque both when using a contra-angle handpiece and when using a ratchet.

Warning

When using surgical instruments with a limited number of permitted uses, this limit must not be exceeded, as doing so may seriously compromise the osseointegration process and, consequently, the treatment outcome.

The maximum insertion torque for dental implants is 50 Ncm. Exceeding the specified maximum insertion torque may cause serious damage to the dental implant and its connection. Refer to the surgical protocol for specific considerations according to implant connection type and bone quality.

Care should be taken not to spill the saline solution when opening the implant vial in those implant models that contain it.

10. Warnings

Treatment planning and the placement of Ziacom® dental implants require specific dental training. Product specifications alone do not guarantee proper use. Users are advised to undertake specialised theoretical and practical training courses in order to become familiar with the techniques and working protocols appropriate for each product, including the biomechanical, radiographic and prosthodontic requirements associated with implant treatment. It is the user's responsibility to remain informed of the current state of the art for each product and its possible applications.

Before using Ziacom® products, it is necessary to be familiar with the corresponding surgical and prosthodontic procedures. The patient must present adequate anatomical and psychological conditions to undergo dental implant treatment. The user must ensure proper treatment planning, guaranteeing treatment safety with a minimal margin of error and respecting vital oral structures and the patient's overall health.

The procedures for the use of surgical drills, taps and other surgical instruments required for implant placement are specified in the surgical protocols corresponding to each implant type. Implant placement and prosthodontic planning must be adapted to the individual conditions of the patient, particularly with regard to the proper distribution of forces. Passive fit of the prosthodontic restoration, correct occlusal adjustment and the avoidance of excessive lateral forces must be achieved. Treatment involving an insufficient number of implants, inappropriate implant size selection or unsuitable positioning to support and transmit the expected loads may result in mechanical failure of the implant, abutment or prosthetic component, clinical screws due to overload or fatigue, and substantial loss of bone surrounding the implant.

Insufficient alveolar bone volume or quality, the occurrence of infection or disease in general, and changes in patient habits are among the potential causes of osseointegration and treatment failure. Deficiencies in hard or soft tissue may result in poor aesthetic outcomes. Inadequate prosthodontic rehabilitation may lead to failure of the entire restorative treatment.

The reuse of single-use products or of any residual liquid contained within the vial, in implant models that contain such liquid, may lead to deterioration of product characteristics. Product geometry may have been affected by previous use or improper handling. Once the packaging has been opened, sterility is lost and reuse entails a risk of surgical or prosthodontic failure and/or harm to patient health.

Electrosurgery is not indicated due to the conductive properties of dental implants.

11. Adverse Effects or Side Effects

Various adverse or undesirable effects may occur in current implant dentistry. These effects are documented in specialised scientific literature and published in journals and textbooks within the dental field.

The most relevant include:

Postoperative discomfort - Local pain - Local inflammation - Haematomas - Haemorrhage - Local or systemic infections - Difficulty speaking - Bone loss and/or fractures - Implant loss - Damage to adjacent teeth - Fracture of implants and/or prosthodontic components - Instrument fracture - Injury to adjacent anatomical structures (inferior alveolar nerve or maxillary sinus).

Ziacom Implants SLU will update these adverse effects whenever necessary, based on information obtained from clinical studies and post-market surveillance data.

12. Information on Liability, Safety and Warranty

The instructions for use and handling of Ziacom® products are based on published international literature, current clinical standards and our clinical experience with our products, and should therefore be understood as general guidance information. Since the handling and use of Ziacom® products are beyond the control of Ziacom Implants SLU, they are the sole responsibility of the user. Ziacom Implants SLU, its subsidiaries and/or its authorised distributors disclaim all liability, whether express or implied, in whole or in part, for any damage or injury resulting from improper handling of the product or from any circumstance not contemplated in its protocols and manuals for the correct use of its products.

The user must ensure that the Ziacom® product used is suitable for the intended procedure and purpose. Neither these instructions for use nor the product handling or working protocols relieve the user of this obligation. The use, handling and clinical application of Ziacom® products must be carried out by qualified professionals holding the required credentials in accordance with the legislation in force in each country. Any use, handling and/or application, whether wholly or partially, of Ziacom® products by unqualified personnel or personnel lacking the necessary qualifications will automatically void any warranty and may cause serious harm to patient health.

Ziacom® products form part of a proprietary system with specific design features and working protocols, including dental implants, prosthetic abutments or components, and surgical or prosthetic instruments. The use of Ziacom® products in combination with elements or components from other manufacturers may result in treatment failure, tissue damage, damage to bone structures, unsatisfactory aesthetic outcomes and serious harm to patient health. For this reason, only original Ziacom® products should be used.

The clinical professional responsible for the treatment is solely responsible for ensuring the use of original Ziacom® products and their use in accordance with the corresponding instructions for use and handling protocols throughout the implant treatment process.

The use of non-original Ziacom® components, instruments or any other products, whether used alone or in combination with any original Ziacom® products, will automatically void the warranty of the original Ziacom® products. Please refer to the Ziacom Implants SLU Warranty Programme (available on the website or by contacting Ziacom Implants SLU, its subsidiaries or authorised distributors).

Any problem or incident related to the device must be reported to the manufacturer, Ziacom Implants SLU, to whom the affected product must also be sent. In the event of a serious incident, the user must additionally submit a report to Ziacom Implants SLU and to the relevant competent authority, as required by local regulations.

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.

Warning

Not all Ziacom® products are available in all countries. Please check availability.

The Ziacom® trademark and other product or service names, as well as their logos, mentioned in this documentation or on the website www.ziacom.com, are registered trademarks of Ziacom Implants SLU.

Ziacom Implants SLU reserves the right to modify, change, discontinue or further develop any products, prices or technical specifications referenced in this document or in any of its catalogues without prior notice. All rights reserved. The total or partial reproduction or publication of this documentation, in any medium or format, is prohibited without the prior written authorisation of Ziacom Implants SLU.

13. Magnetic Resonance Imaging (MRI)

The safety and compatibility of Ziacom® products in a magnetic resonance (MR) environment have not been evaluated. They have not been tested for heating or migration in such an environment.

14. Disposal of Materials

Dental implants, abutments and instruments, as well as any product and consumable used during dental implant placement surgery and subsequently discarded, may pose a risk to the health of individuals handling them. Therefore, it is recommended to consult and comply with the applicable legislation on biohazardous waste.

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Titansure

EN Carefully remove the seal from the implant blister.

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
Titansure
Active

EN Hold the container in a horizontal position and remove the cap by turning it in the opposite direction.

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Titansure


EN Let the vial carrying the implant fall onto a sterile drape in the surgical field.

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
Titansure

EN Hold the vial with one hand in a vertical position. Remove the cap by turning it vertically.

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
Titansure

EN Hold the implant carrier vial with one hand and put in the insert, for ratchet or CA, into the ZPlus Mount.

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
Titansure

EN Remove the implant–plastic carrier assembly by pulling it vertically upwards from the vial.

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Titansure
Active

EN Hold the vial with one hand in a vertical position. Remove the cap by turning it vertically.
Note: Exercise caution when opening the vial, as the implant is immersed in a saline solution.

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


EN

Hold the implant vial vertically with one hand and insert the manual or handpiece adapter into the implant:
Option 1 – With ratchet / Option 2 – With contra-angle handpiece / Option 3 – With manual adapter

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
EN

Insert the iDrive mechanical adapter into the implant head, aligning the hexagon of the carrier with the internal hexagon of the implant, and press until fully seated. Only for ZM8 N and ZM8 S No Mount implants.

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
EN

Remove the implant by pulling it vertically upwards from the vial:
Option 1 – With ratchet / Option 2 – With contra-angle handpiece / Option 3 – With manual adapter

Note: Do not reuse any remaining saline solution.

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Ziacom Medical USA LLC