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EN | INSTRUCTIONS FOR USE

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IMPLANT INSTRUCTIONS
FOR USE

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Ziacom® orthodontic implants

INSTRUCTIONS FOR USE

Instructions for use applicable to the DSQ orthodontic implants systems.

IMPORTANT INFORMATION

Read these instructions for use carefully before using the product.

This document contains basic information for the use of original Ziacom® orthodontic implants, micro-implants or mini-screws, hereinafter "Ziacom® orthodontic implants" or simply "Ziacom® products". This document has been written as a quick reference guide for the physician in charge of the treatment of orthodontia, henceforth, "User". It does not provide sufficient indications and technical specifications for the use of Ziacom® products. It is neither an alternative nor a substitute for specialised training and professional clinical experience.

Ziacom® products must be used with proper treatment planning and following the working protocols established by Ziacom Medical SL. Read carefully: the specific working protocols for this product, as well as the instructions for use and the instructions for cleaning, disinfection and sterilisation before using this Ziacom® product. You can consult them on our website www.ziacom.com or request them from your official distributor authorised by Ziacom Medical SL.

1. Description of the system

Ziacom® products consist of orthodontic implants, dental implants, prosthetic abutments or components and surgical and prosthetic instruments.

Ziacom Medical SL has created two different types of orthodontic implants to serve as temporary anchor points in orthodontic treatment according to the existing scientific literature and current clinical standards. Ziacom Medical SL develops, manufactures and commercialises latest-generation implant products and a wide range of restorative solutions to respond to the particular needs of each patient.

2. Description of the products

Orthodontic implants

Ziacom® orthodontic implants are manufactured in biocompatible grade 5 titanium alloy, in accordance with ASTM F136-13 and UNE-EN ISO 5832-3:2022 standards, without surface treatment to avoid excessive osseointegration that hinders explantation at the end of treatment. Ziacom® orthodontic implants are manufactured in different diameters, lengths, heads (cross or button type) and collars depending on their intended purpose or location.

3. Indications for use

It is essential to previously evaluate the patient by means of a clinical and imaging diagnosis and study of models, as well as to determine the possible contraindications for the treatment, prior to carrying out a treatment with Ziacom® products. A correct orthodontic treatment that requires the use of orthodontic implants must start from a diagnosis, planning and design of orthodontic treatment that meets the needs and expectations of the patient.

Ziacom® orthodontic screws serve as a fixation point, which operates as an anchorage point to facilitate the various orthodontic treatments in which it is used. They are supplied with different head designs and in different lengths and diameters so that they can be adapted to different clinical cases.

Ziacom® implants are single-use products. The head design of orthodontic implants allows their use together with other products and accessories commonly used in orthodontics (orthodontic wires, steel ties, elastic ligatures, brace chains, springs, etc.), for more information consult the catalogue.

The insertion and removal of the orthodontic implant must be performed by a qualified professional using the recommended instruments and following the indications in the catalogues and working protocols.

Carefully consult the clinical indications for each orthodontic implant model in the relevant catalogues and working protocols before planning treatment. Do not use Ziacom® orthodontic implants in locations or for purposes that are not described as indications in your catalogues or working protocols.

4. Contraindications

A preoperative medical examination of the patient is necessary to determine whether there are absolute and/or relative contraindications or risk factors for treatment with orthodontic implants. It is the user's responsibility to evaluate the potential benefits and risks of using orthodontic implants in orthodontic treatment and to make appropriate decisions on a case-by-case basis. Ziacom® products should not be used in patients who lack the minimum medical conditions for orthodontic treatment and restoration with orthodontic implants.

General contraindications:

Patients with systemic diseases (disorders in blood coagulation, uncontrolled endocrine disorders, rheumatic diseases, liver diseases, neoplasms, diabetes or other serious diseases), psychological or psychic alterations, poor oral hygiene, severe or uncontrolled periodontal disease, inadequate quantity/quality of bone that prevents the retention of the orthodontic implant in the bone cortex, inadequate oral habits (aggressive tooth brushing, smoking, foods with high sugar content, oral piercings...).

Ziacom Medical SL will update these contraindications when necessary, based on the information obtained through clinical studies and post-sales information.

5. Intended patients

The target population for these medical devices are those patients who, in the opinion of the clinician assessing their case, require treatment of the mandible and maxilla, and are intended to restore chewing function, as well as provide an improvement in occlusion and the distribution of occlusal loads, phonation, aesthetics, prosthetic stability and maintenance of the alveolar process, as well as a positive psychological effect.

The recommended age to initiate an orthodontic treatment is 6 years old.

6. Storage

Ziacom® products must be stored in their original packaging, in a clean, dry place, protected from sunlight and adverse conditions. The storage temperature must be lower than +55°C and higher than -10°C.

7. Packaging and supply state

Ziacom® orthodontic implants are supplied in a cardboard box with a product identification label, which acts as a security seal against handling.

This package contains:

- This explanatory product information leaflet in paper or electronic format.
- Thermo-sealed plastic blister with removable labels on the back, indicating the manufacturer, complete product reference, lot number, dimensions of the orthodontic implant and the relevant symbol (explained in these Instructions for use). Inside is the orthodontic implant.

The removable labels serve as a traceability control of the implant. This is achieved by sticking one label to the implant card for the patient, and the rest for the clinic's record.

Ziacom Medical SL guarantees that all its products follow an exhaustive process of manufacture, control and cleaning before being packaged. Before using Ziacom® products you must check the correct integrity of the package and make sure that it is not damaged. If there is any defect or damage, it must not be used and must be notified to Ziacom Medical SL or to the official authorised distributor.

8. Cleaning and sterilisation

Ziacom® orthodontic implants are supplied non-sterile.

Orthodontic implants that will come into contact with the mouth must be cleaned, disinfected and sterilised prior to use following the cleaning, disinfection and sterilisation protocol for orthodontic implants, instruments, surgical boxes and abutments of Ziacom Medical SL. You can consult it in our web page www.ziacom.com or in the catalogues of Ziacom Medical SL.

9. Precautions

Insertion and instruments

The working protocols describe in detail the important precautions and factors to be taken into account during the orthodontic implant placement process. Before insertion, it is essential to consult the working protocol of the product, as well as all indications for its correct use in the Ziacom® product catalogues.

Special attention must be paid and the necessary measures must be taken to prevent the orthodontic implant or instruments from being swallowed or bronchoaspirated by the patient.

Pay special attention not to exert excessive forces or torques during the insertion of the orthodontic implant. If during insertion you detect a variation in bone consistency, stop insertion immediately and make sure that you are not inserting into inappropriate bone or dental structures such as dental roots.

Attention

The maximum insertion depth of orthodontic implants is up to the collar of the orthodontic implant. Avoid applying excessive insertion torque to orthodontic implants as this can cause serious damage to the orthodontic implant and bone. The application of the insertion torque should be under the responsibility of the user.

10. Warnings

Treatment planning and placement of Ziacom® orthodontic implants require specific dental training. Product specifications alone do not guarantee proper use. Users are encouraged to take specialised theoretical and practical training courses to understand the appropriate working techniques and protocols for each product, including biomechanical, radiographic and prosthetic requirements associated with orthodontic treatments.

It is the user's obligation to know and be informed of the current state of the technique for each product and its possible applications.

Before using Ziacom® products it is necessary to know and be familiar with the corresponding working procedures. The patient must have sufficient anatomical and psychological conditions to undergo dental treatment. The user must ensure the correct planning of the treatment, guaranteeing its safety with a minimum margin of error and respecting the vital structures of the mouth and the general health of the patient.

The procedure for the use of the insertion instruments required for the placement of orthodontic implants is specified in the working protocols for each type of orthodontic implant. Orthodontic implant placement and orthodontic planning should be tailored to individual patient conditions. The point of insertion should be chosen so that the load should have a 90° angle or be perpendicular to the longitudinal axis of the orthodontic implant. Under or overloaded treatment, an inappropriate choice of size or an inappropriate position to support the expected loads can lead to failure of the orthodontic implant or treatment (overload or fatigue).

A qualitatively or quantitatively insufficient alveolar bone, the appearance of infections or diseases in general, and changes in the habits of the patient, are some potential causes of failure of the treatment. Lack of bone or soft tissue can lead to poor outcome and inadequate loading can lead to failure of the entire restorative treatment.

The reuse of single-use products carries the danger of possible deterioration of their properties. The geometry of the product may have suffered damage during previous use or due to an unsuitable use. This implies risk of orthodontic failure and/or damage to the patient's health.

Electrosurgery is not indicated due to the conductivity of orthodontic implants.

11. Adverse or side effects

Many adverse or undesired effects can present themselves in orthodontic. These are documented in the specialised scientific bibliography and published in journals and books specialised in the dentistry sector.

However, the most relevant are:

Operative or postoperative pain, inflammation/infection of the mucosa of the area, haematomas or haemorrhages, bone mobility, loss or fall of the orthodontic implant, alveolar fracture, fracture of the orthodontic implant, damage to adjacent roots or teeth, injury to adjacent anatomical structures (nostrils, maxillary sinus, vascular/nerve structures).

Ziacom Medical SL will update these side effects when necessary, based on the information obtained through clinical studies and post-sales information.

12. Information about responsibility, safety and guarantee

The handling and use indication of Ziacom® products are based on published international literature, current clinical standards and our clinical experience with our products and should therefore be understood as indicative general information. The handling and use of Ziacom® products, as they are beyond the control of Ziacom Medical SL, are the sole responsibility of the user. Ziacom Medical SL, its subsidiaries and/or its authorised official distributors decline all responsibility, express or implicit, total or partial, for the possible damages caused by the bad handling of the product or by any other fact not contemplated in its protocols and instructions for the correct use of its products.

The user of the product must ensure that the Ziacom® product used is suitable for the procedure and the intended purpose. Neither these instructions for use, nor the working protocols of the products, exempt the user from this obligation. The use, handling and clinical application of Ziacom® products must be carried out by qualified professional personnel with the necessary certification according to the legislation in force in each country. The use, manipulation and/or application, in whole or in part, of Ziacom® products in any of their manufacturing phases by unqualified personnel or without the necessary qualifications, automatically voids any type of guarantee and may cause severe damage to the health of the patient.

Ziacom® products have their own design features and working protocols, so the insertion of orthodontic implants must be done with the adapters and insertion instruments recommended by Ziacom Medical SL in the catalogue and in the working protocols. However, the design of orthodontic implants enables them to be used together with other products and accessories commonly used in orthodontics (orthodontic wires, steel ties, elastic ligatures, brace chains, springs, etc.). Insertion of Ziacom® orthodontic implants with adapters or instruments from other manufacturers may result in insertion failure, deterioration of the orthodontic implant, injury to the insertion site by causing tissue damage, bone structure damage, inadequate aesthetic results and serious damage to the patient's health. For this reason, insertion should only be performed with instruments and adapters recommended by Ziacom Medical SL.

The clinician in charge of treatment is the sole responsible for ensuring that the original Ziacom® products are used and that they are used in accordance with the corresponding instructions for use and working protocols throughout the orthodontic treatment process. The use of non-original Ziacom® instruments used alone or in combination with any of the original Ziacom® products will automatically void the guarantee of the original Ziacom® products. Consult the Ziacom Medical SL Guarantee Programme (available on the website or by contacting Ziacom Medical SL, its subsidiaries or authorised distributors).

Any problem or incident with regard to the device should be reported to the manufacturer, Ziacom Medical S.L., to which the user should also send the affected product. In the event of a serious incident, the user should also submit a report to the manufacturer, Ziacom Medical S.L. as well as the relevant competent authority, according to the requirements of local legislation.

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.

Warning

Not all Ziacom® products are for sale in every country. Check availability in your country.

The Ziacom® brand and the names of other products and services, including their logos, that are mentioned in this document or on the website www.ziacom.com, are registered trademarks of Ziacom Medical S.L.

Ziacom Medical S.L. reserves the right to modify, change, eliminate or update any of the products, prices or technical specifications mentioned in this document or in any of its catalogues without prior notification. All rights reserved. The reproduction of this document, whole or in part and in any medium or format, without the corresponding written authorisation from Ziacom Medical SL is prohibited.

13. Magnetic Resonance

The safety and compatibility of Ziacom® products in a magnetic resonance environment has not been evaluated, nor have they undergone testing to verify whether they heat up or migrate in such environments.

14. Removal of the material

The dental implants, abutments and kit, as well as each product and consumable used during surgery for placing dental implants and which are disposed of may pose a health hazard to people handling them. For this reason, you are advised to check and comply with current legislation on biohazardous waste.

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EN | CLEANING, DISINFECTION AND STERILISATION

Cleaning, disinfection and sterilisation

The protocols described in this section must only be carried out by personnel qualified to clean, disinfect and sterilise the dental materials specified here in.

CLEANING AND DISINFECTION INSTRUCTIONS

Applicable for instruments, surgical and prosthetic boxes.

Disassembly

1. Dismount* the appropriate instruments, for example manual ratchets, drills or drill stops.
2. Remove the various components from the surgical or prosthetic box for correct cleaning.

Cleaning and disinfection

For disinfecting instruments and surgical boxes:

1. Submerge the instruments in a detergent/disinfectant solution** suitable for dental instruments to help eliminate any adhered biological residues. If an ultrasound bath is available***, confirm that the detergent/disinfectant solution is indicated for use with this type of equipment.
2. Manually remove any biological residues with a non-metallic brush and pH-neutral detergent.
3. Rinse with copious water.
4. When cleaning the surgical and prosthetic boxes, always use a pH-neutral detergent and non-abrasive utensils to avoid damaging the surface of the boxes.
5. Dry the materials with disposable cellulose, lint-free clothes or compressed air.

For disinfecting plastic caps and spacers:

1. Submerge in a neat benzalkonium chloride solution for 10 minutes.
2. Rinse with distilled water.
3. Dry the caps and spacer before use.

Inspection

1. Check that the instruments are perfectly clean; if not, repeat the cleaning and disinfection steps.
 2. Discard any instruments with imperfections and replace them before the next procedure.
 3. Check that the instruments and the surgical and prosthetic boxes are perfectly dry before reassembling the parts and proceeding to their sterilisation.
- * See assembly and disassembly instructions.
- ** Follow the instructions from the disinfectant's manufacturer to determine the correct concentrations and times.
- *** Follow the instructions from the ultrasound bath's manufacturer to determine the correct temperature, concentration and times.

STERILISATION INSTRUCTIONS FOR STEAM AUTOCLAVE

Applicable to orthodontic implants, abutments, instruments, surgical and prosthetic boxes, pins, fixation screws and laminates.

1. Introduce each material separately in individual sterilisation bags, then seal the bags. For joint sterilisation, place the instruments in their surgical box, introduce the box into a sterilisation bag and seal the bag.
2. Place the bags to be sterilised in the autoclave.
3. Sterilise in a steam autoclave at 134°C/273°F (max. 137°C/276°F) for 4 min (minimum) and at 2 atm. Torque wrenches must be sterilised in 3 vacuum cycles at 132°C/270°F for a minimum of ≥ 4 minutes and vacuum-dried for a minimum of 20 minutes.

For the United States only: The validated and recommended sterilisation cycle for the US must be performed in a steam autoclave at 132°C/270°F for at least 15 min and with the drying time of at least 15- 30 min.

Important:

- Make sure the drying stage is allowed to run to completion, otherwise the products may be damp.
- Check the sterilisation equipment if the materials or sterilisation bags are damp at the end of the sterilisation cycle.
- Perform the necessary maintenance actions on the autoclave according to the established periodicity and following the manufacturer's instructions.

STORAGE OF ZIACOM® PRODUCTS

- Store the products in their original packaging and in a clean, dry location until they are used.
- After sterilisation, keep the products in the sealed sterilisation bags and in a clean, dry location.
- Never exceed the use by date indicated by the manufacturer of the sterilisation bags.
- Always follow the indications of the manufacturer of the sterilisation bags.

GENERAL RECOMMENDATIONS

- Never use damaged or dirty material; never reuse single-use products. The user is responsible for following the instructions described in this document correctly.
- The attention to piercing or sharp elements. Gloves should be worn when cleaning the materials to avoid accidents during handling.
- Follow the safety instructions indicated by the manufacturer of the disinfectant agent.
- The product's sterility cannot be guaranteed if the sterilisation bag is open, damaged or damp.
- Respect all stages of the sterilisation process. If the materials or sterilisation bags contain traces of water or moisture, check the autoclave and repeat the sterilisation.
- Orthodontic abutments and implants are supplied UNSTERILISED and must always be sterilised before use.
- Instruments and surgical and prosthetic boxes are supplied UNSTERILISED and must always be sterilised before use and cleaned and disinfected after use.
- The sterilisation, cleaning and disinfection processes gradually deteriorate the instruments. Inspect the instruments thoroughly to detect any signs of deterioration.
- Avoid contact between products made from different materials (steel, titanium, etc.) during the cleaning, disinfection and sterilisation processes.
- Ziacom Medical SL recommends these instructions are implemented for the correct maintenance and safety of their products; accordingly, the company refuses any liability for any damage to the products that could arise if the user applies alternative cleaning, disinfection and sterilisation procedures.

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EN Press on the word “PRESS” and tear open the cardboard box. Remove the cardboard flap and take out the blister.

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EN Carefully remove the seal from the blister.

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EN Remove the plastic holder containing the DSQ orthodontic implant from the blister pack.

Important note: DSQ implants must not be sterilised in their original packaging or with the plastic holder.

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
EN Remove the DSQ orthodontic implant from its plastic holder by taking it by the transmucosal collar using pliers and place it in the DSQ implant surgical box. Sterilise this box according to the instructions for use.

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EN Once the instruments have been sterilised, use the insertion key in the handle or the counter-angle insertion key to hold the implant for clinical use.

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EN Remember to remove the identification labels from the implant to stick them on the Implant Card and on the patient history sheet in order to document the traceability of the reference and batch of the implant.

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Ziacom Medical USA, LLC
333 S.E. 2nd Avenue, Suite 2000
Tel: +1 (786) 224-0089 - Miami, FL 33131 - USA