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EN | INSTRUCTIONS FOR USE

IMPLANT INSTRUCTIONS
LAMINATES

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Ziacom® ZS3 Laminates

INSTRUCTIONS FOR USE

Instructions for use laminates Ziasure ZS3.

IMPORTANT INFORMATION

Please read these instructions carefully before using the product.

This document contains essential information on the use of Ziacom® original osteosynthesis screws, laminates, hereinafter "Ziacom® laminates" or simply "Ziacom® products". This documentation has been prepared as a quick guide to be consulted by the clinician responsible guided bone regeneration, hereinafter the "user", and does not provide sufficient guidelines or technical specifications for correct use of ZIACOM® products. This is not an alternative nor substitute for specialised training or professional clinical experience.

Ziacom® products must be used according to a suitable treatment plan and in strict compliance with the working protocols established by Ziacom Medical SL Please read the specific working protocols for this product, as well as the instructions for use and cleaning, disinfection and sterilisation carefully before using this Ziacom® product. You can find this information on our website at www.ziacom.com or request it from your Ziacom Medical SL authorised distributor.

1. Description of the System

This Ziacom® product is composed of fixation laminates.

Ziacom Medical SL has created the ZS3 Ziasure titanium laminates whose design acts as physical barriers in guided bone regeneration procedures, creating and maintaining the necessary space for bone formation and providing stability to the coagulum and bone substitute, which translates into predictable guided bone regeneration procedures. Ziacom Medical SL also develops, manufactures and markets state-of-the-art implantological products and a wide range of restorative solutions to meet the particular needs of each patient.

2. Description of the Products

Ziacom® fixation laminates are made of biocompatible grade 2 titanium alloy, in accordance with ASTM F67 and UNE-EN ISO 5832-2:2018 standards, without surface treatment to avoid excessive osseointegration that hinders explantation at the end of treatment.

3. Indications for Use

Ziacom® laminates are indicated for the stabilisation and temporary fixation of membranes, biomaterials used in guided bone regeneration (GBR) or orthodontic appliances in the maxilla or mandible.

Before performing a treatment with Ziacom® products, it is essential to previously evaluate the patient by means of a clinical and imaging diagnosis and a study of models, as well as to determine the possible contraindications for such treatment. A correct bone regeneration treatment facilitates the formation of viable bone tissue for the placement of a dental implant or to treat some type of bone defect.

4. Contraindications

The patient must undergo a preoperative medical examination, as well as assessment of the surgical conditions of the area and of the procedure to determine whether there are any absolute and/or relative contraindications or risk factors for the use of fixation pins. The user is responsible for evaluating the potential risks and benefits of the use of the fixation pins in bone regeneration procedures. Ziacom® products should not be used in patients who lack the minimal medical conditions necessary to carry out regeneration treatment. Neither should dental implants and prosthetic rehabilitations using dental implants be placed in patients who do not meet the minimum medical conditions.

General Contraindications:

Patients with systemic diseases (blood clotting disorders, endocrine disorders, rheumatic diseases, liver disease, neoplasms or other serious diseases), psychological or mental disorders, poor oral hygiene, untreated severe periodontal disease, toxic habits (smoking, drug addiction, alcoholism), active infection or pathological processes at the site of placement, allergies to the fixation pin materials, anatomical structures at a distance of less than 2 mm from the site of insertion (teeth, inferior alveolar canal, nasopalatine duct, maxillary sinus, etc.)

The user should review the compatibility of certain medications with bone regeneration treatment.

Ziacom Medical SL will update these contraindications as necessary, based on information obtained from clinical studies and/or post-marketing information.

5. Intended patients

The target population for these medical devices are those patients who, in the opinion of the clinician assessing their case, require treatment in the mandible and maxilla, and are intended to restore chewing function, as well as providing an improvement in occlusion and the distribution of occlusal loads, phonation, aesthetics, prosthetic stability and maintenance of the alveolar process, as well as a positive psychological effect.

These patients must have completed their growth. The age guideline is older than 18 years, this can be confirmed after cessation of growth with a wrist X-ray. Patients at an early age patients will be specifically informed of the changes they may need to make to the prostheses during their lifetime with a with the consequent aesthetic and economic impact due to functional ankylosis of the implant in living tissue.

6. Storage

Ziacom® products must be stored in their original packaging in a dry, clean place away from sunlight and adverse conditions. They must be stored at a temperature below +55°C and above -10°C.

7. Packaging and Supply Format

Ziacom® laminates are supplied in a cardboard box with a product identification label, which functions as a security seal to protect against tampering.

The pack contains:

- The QR code to access these explanatory instructions for use of the product.
- Plastic heat-sealed blister and removable labels on the back, indicating the manufacturer, full product reference, batch number, laminates dimensions and relevant symbols (explained in these instructions for use). The laminates is enclosed inside.

The removable labels are intended for traceability of the fixation screws, by sticking one label on the implant card of the patient and the remaining ones on the clinical record.

Ziacom Medical SL guarantees that all of its products follow a strict manufacturing, monitoring and cleaning process before being packaged. Before using Ziacom® products, you must check that the packaging is intact and ensure that it is not damaged. If there is any damage or defect, the product must not be used and Ziacom Medical SL or your authorised official distributor must be notified.

8. Cleaning and Sterilisation

Ziacom® laminates are supplied NON sterile.

Laminates that are to come into contact with the mouth must be cleaned, disinfected and sterilised before use, following the Ziacom Medical SL cleaning, disinfection and sterilisation instructions. You can find this on our website at www.ziacom.com or in Ziacom Medical SL catalogues.

9. Precautions

Insertion and Instruments

The working protocols describe in detail the precautions and important factors to be taken into consideration during placement of the laminates. Prior to placement, it is essential to have consulted the working protocol for the product, along with all other indications regarding the correct use as indicated in the Ziacom® product catalogues.

Special care must be taken to take the necessary steps to prevent the laminates or instruments from being swallowed or aspirated by the patient.

Drilling of the graft in place must be performed while always controlling the drilling axis, with tapping, irrigation and a slow revolution speed, penetrating a maximum of 3 mm of cortical bone in the recipient bed.

Take extra care not to apply excessive force during placement of the laminates. If during placement a variation is detected in the bone consistency, stop insertion immediately and make sure that you are not inserting into inadequate bone or dental structure, for example into tooth roots, as well as if you detect any deformity in the body of the laminates that may lead to any type of fracture. It is recommended that a maximum of 5 Ncm torque in 1.00 mm and 1.20 mm diameter screws, and 15 Ncm in 1.50 mm and 2.00 mm diameter screws is not exceeded.

Please note

The correct insertion depth of the laminates is achieved when the graft, material and/or membrane is secured and stable, making sure that it is immobile.

10. Warnings

Planning regenerative treatment, selection of the biomaterials to use and placement of Ziacom® laminates requires specific training in dentistry. The product specifications themselves are no guarantee of proper use. We recommend that users take specialised theoretical and practical training courses in order to understand the proper working techniques and protocols for each product, including the biomechanical, imaging and biological requirements associated with bone regeneration treatment.

It is the user's responsibility to know and find out about the current status of the technique for each product and its possible applications.

It is necessary to be familiar with and understand the corresponding working procedures before using Ziacom® products. The patient must meet sufficient medical, anatomical and psychological conditions to undergo treatment with dental implants. The user must ensure proper treatment planning is carried out, to guarantee the safety of the treatment with a minimal margin of error, respecting the patient's vital mouth structures and general health.

The procedure for using the insertion kit for placing the laminates is specified in the corresponding working protocols for each type of laminates. Placement of the laminates and the bone regeneration procedures must be tailored to the conditions of the individual patient.

Ziacom® laminates are single-use products. Reuse of single-use products carries the risk of possible deterioration of their characteristics. The geometry of the product may have suffered damage from previous use or due to inadequate use. This entails the risk of failure of the bone regeneration and/or damage to the patient's health.

11. Adverse Events or Side Effects

Regenerative treatments may lead to various adverse or undesirable effects that are documented in the specialised scientific literature published in specialised journals or books on dentistry.

However, the most relevant are:

Laceration of the membrane, injury to adjacent anatomical structures (nasal passages, maxillary sinus, vascular/nerve structures), infections, entrapment of the screw, etc.

Ziacom Medical SL will update these side effects as necessary, based on information obtained from clinical studies and post-marketing information.

12. Liability, safety and warranty information

The instructions for use and handling of Ziacom® products are based on internationally published literature, current clinical standards and our own clinical experience so should be understood as general guiding information. The handling and use of Ziacom® products is the sole responsibility of the user as it is outside the control of Ziacom Medical SL. Ziacom Medical SL, their subsidiaries and/or their official authorised distributors do not accept any liability, whether explicit or implicit, total or partial, for possible damage or injury caused by improper handling of the product or any other situation not considered in their protocols and instructions for the correct use of their products.

The user must ensure that the Ziacom® product is appropriate for the intended procedure and end purpose. Neither these instructions for use nor the working or handling protocols for the products release the user from this obligation. Ziacom® products must be used, handled and applied clinically by professionals with the appropriate training and qualifications required according to current legislation in each country. Total or partial use, handling or application of any phase of Ziacom® products by unqualified personnel, or by personnel that does not have the necessary training, will automatically void any type of warranty and may cause serious damage to the patient's health.

Ziacom® products have their own design features and working protocols, so the insertion of the laminates should be carried out with the adaptors and insertion kit recommended by Ziacom Medical SL in the catalogue and the working protocols. Insertion of Ziacom® laminates with adaptors or instruments from other manufacturers may lead to insertion failure, deterioration of the laminates, adverse impact on the area of insertion in causing damage to the tissues and bone structures. For this reason, insertion should be carried out with the kit and adaptors recommended by Ziacom Medical, SL

The design of the laminates enables them to be used together with other products, such as various types of membranes and biomaterials.

The clinician in charge of the treatment is solely responsible for ensuring that original Ziacom® products are used and that they are used according to the corresponding instructions for use and working protocols throughout the entire

guided bone regeneration treatment. The use of any other non-original Ziacom® instruments, whether alone or in combination with any original Ziacom® products, will immediately void the warranty of the original Ziacom® products. See the Ziacom Medical SL Warranty Programme (available on the website or by contacting Ziacom Medical SL, their subsidiaries or authorised distributors).

Any problem or incident with regard to the device should be reported to the manufacturer, Ziacom Medical SL, to which the user should also send the affected product. In the event of a serious incident, the user should also submit a report to the manufacturer, Ziacom Medical SL as well as the relevant competent authority, according to the requirements of local legislation.

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.

Warning

Not all Ziacom® products are for sale in every country. Check availability in your country.

The Ziacom® brand and the names of other products and services, including their logos, that are mentioned in this document or on the website www.ziacom.com, are registered trademarks of Ziacom Medical SL

Ziacom Medical SL reserves the right to modify, change, eliminate or update any of the products, prices or technical specifications mentioned in this document or in any of its catalogues without prior notification. All rights reserved. The reproduction of this document, whole or in part and in any medium or format, without the corresponding written authorisation from Ziacom Medical SL is prohibited.

13. Magnetic Resonance

The safety and compatibility of Ziacom® products in a magnetic resonance environment has not been evaluated, nor have they undergone testing to verify whether they heat up or migrate in such environments.

14. Removal of the material

The dental implants, abutments and kit, as well as each product and consumable used during surgery for placing dental implants and which are disposed of may pose a health hazard to people handling them. For this reason, you are advised to check and comply with current legislation on biohazardous waste.

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EN | CLEANING, DISINFECTION AND STERILISATION

Cleaning, disinfection and sterilisation

The protocols described in this section must only be carried out by personnel qualified to clean, disinfect and sterilise the dental materials specified here in.

CLEANING AND DISINFECTION INSTRUCTIONS

Applicable for instruments, surgical and prosthetic boxes.

Disassembly

1. Dismount* the appropriate instruments, for example manual ratchets, drills or drill stops.
2. Remove the various components from the surgical or prosthetic box for correct cleaning.

Cleaning and disinfection

For disinfecting instruments and surgical boxes:

1. Submerge the instruments in a detergent/disinfectant solution** suitable for dental instruments to help eliminate any adhered biological residues. If an ultrasound bath is available***, confirm that the detergent/disinfectant solution is indicated for use with this type of equipment.
2. Manually remove any biological residues with a non-metallic brush and pH-neutral detergent.
3. Rinse with copious water.
4. When cleaning the surgical and prosthetic boxes, always use a pH-neutral detergent and non-abrasive utensils to avoid damaging the surface of the boxes.
5. Dry the materials with disposable cellulose, lint-free clothes or compressed air.

For disinfecting plastic caps and spacers:

1. Submerge in a neat benzalkonium chloride solution for 10 minutes.
2. Rinse with distilled water.
3. Dry the caps and spacer before use.

Inspection

1. Check that the instruments are perfectly clean; if not, repeat the cleaning and disinfection steps.
 2. Discard any instruments with imperfections and replace them before the next procedure.
 3. Check that the instruments and the surgical and prosthetic boxes are perfectly dry before reassembling the parts and proceeding to their sterilisation.
- * See assembly and disassembly instructions.
- ** Follow the instructions from the disinfectant's manufacturer to determine the correct concentrations and times.
- *** Follow the instructions from the ultrasound bath's manufacturer to determine the correct temperature, concentration and times.

STERILISATION INSTRUCTIONS FOR STEAM AUTOCLAVE

Applicable to orthodontic implants, abutments, instruments, surgical and prosthetic boxes, pins, fixation screws and laminates.

1. Introduce each material separately in individual sterilisation bags, then seal the bags. For joint sterilisation, place the instruments in their surgical box, introduce the box into a sterilisation bag and seal the bag.
2. Place the bags to be sterilised in the autoclave.
3. Sterilise in a steam autoclave at 134°C/273°F (max. 137°C/276°F) for 4 min (minimum) and at 2 atm. Torque wrenches must be sterilised in 3 vacuum cycles at 132°C/270°F for a minimum of ≥ 4 minutes and vacuum-dried for a minimum of 20 minutes.

For the United States only: The validated and recommended sterilisation cycle for the US must be performed in a steam autoclave at 132°C/270°F for at least 15 min and with the drying time of at least 15- 30 min.

Important:

- Make sure the drying stage is allowed to run to completion, otherwise the products may be damp.
- Check the sterilisation equipment if the materials or sterilisation bags are damp at the end of the sterilisation cycle.
- Perform the necessary maintenance actions on the autoclave according to the established periodicity and following the manufacturer's instructions.

STORAGE OF ZIACOM® PRODUCTS

- Store the products in their original packaging and in a clean, dry location until they are used.
- After sterilisation, keep the products in the sealed sterilisation bags and in a clean, dry location.
- Never exceed the use by date indicated by the manufacturer of the sterilisation bags.
- Always follow the indications of the manufacturer of the sterilisation bags.

GENERAL RECOMMENDATIONS

- Never use damaged or dirty material; never reuse single-use products. The user is responsible for following the instructions described in this document correctly.
- The attention to piercing or sharp elements. Gloves should be worn when cleaning the materials to avoid accidents during handling.
- Follow the safety instructions indicated by the manufacturer of the disinfectant agent.
- The product's sterility cannot be guaranteed if the sterilisation bag is open, damaged or damp.
- Respect all stages of the sterilisation process. If the materials or sterilisation bags contain traces of water or moisture, check the autoclave and repeat the sterilisation.
- Orthodontic abutments and implants are supplied UNSTERILISED and must always be sterilised before use.
- Instruments and surgical and prosthetic boxes are supplied UNSTERILISED and must always be sterilised before use and cleaned and disinfected after use.
- The sterilisation, cleaning and disinfection processes gradually deteriorate the instruments. Inspect the instruments thoroughly to detect any signs of deterioration.
- Avoid contact between products made from different materials (steel, titanium, etc.) during the cleaning, disinfection and sterilisation processes.
- Ziacom Medical SL recommends these instructions are implemented for the correct maintenance and safety of their products; accordingly, the company refuses any liability for any damage to the products that could arise if the user applies alternative cleaning, disinfection and sterilisation procedures.

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EN Press the word "PRESS" and open the carton. Remove the flap from the carton and take out the accessory bag.

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EN Open the accessory bag at the easy-open area to remove the ZS3 laminate.

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

EN Hold the ZS3 laminate with tweezers and place them in the ZS2 surgical kit box. Sterilise the box according to the instructions for use.

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
EN Once the instruments have been sterilised, use a forceps to pick up the ZS3 laminates for clinical use.

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EN Remember to remove the identification labels from the ZS3 laminates to stick them on the Implant Card and on the patient's clinical records, thereby enabling traceability of the product.

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Ziacom Medical USA, LLC
333 S.E. 2nd Avenue, Suite 2000
Tel: +1 (786) 224-0089 - Miami, FL 33131 - USA