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EN | INSTRUCTIONS FOR USE

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FR | INSTRUCTIONS D'USAGE

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Di²gital Arch® Ziacom®

INSTRUCTIONS FOR USE

Instructions for use applicable to the implant systems impression screw Di²gital Arch® 1.0.

IMPORTANT INFORMATION

Please read these instructions carefully before using the product.

This document contains basic information on the use of the Ziacom® original Di²gital Arch® impression screw, hereinafter referred to as the Ziacom Di²gital Arch® system, or simply “Ziacom® products”. This documentation has been prepared as a quick guide to be consulted by the clinician responsible for treatment, hereinafter the “user”, and does not provide sufficient indications or technical specifications for correct use of Ziacom® products. This is not an alternative nor substitute for specialised training or professional clinical experience.

Ziacom® products must be used according to a suitable treatment plan and in strict compliance with the manufacturer’s working protocols: Ziacom Implants SLU Please read the specific working protocols for this product, as well as the instructions for use and cleaning, disinfection and sterilisation carefully before using this Ziacom® product. You can find this information on our website at www.ziacom.com or request it from your Ziacom Implants SLU authorised distributor.

1. Description of the System

Ziacom® products consist of a Di²gital Arch® impression screw, digital impression abutments and surgical kit.

Ziacom Implants SLU has created two Di²gital Arch® 1.0 impression screws which serve as a reference structure to facilitate alignment of virtual files in various clinical situations in which you wish to copy the prior situation of the patient in prosthetic procedures with intraoral scanners, in accordance with the existing scientific literature and current clinical standards. Ziacom Implants SLU also develops, manufactures and markets latest-generation implant surgery products and a wide range of restoration solutions to meet each patients particular needs.

2. Description of the Products

The Ziacom® Di²gital Arch® 1.0 impression screw is manufactured in biocompatible grade-5 titanium alloy, in accordance with standards ASTM F136-13 and UNE-EN-ISO 5832-3:2012, without surface treatment. Ziacom® impression screws are manufactured with the same head in different lengths.

3. Indications for Use

Before carrying out treatment with Ziacom® products, it is vital to assess the patient beforehand with clinical diagnosis, imaging and study models, as well as determine possible contraindications for the treatment. Correct impression taking involves creating a prosthesis and planning the restoration treatment, addressing functional defects (phonation, mastication etc.).

The Ziacom® Di²gital Arch® 1.0 is indicated to serve as a reference structure when aligning the virtual files of a single patient at different stages of treatment or in different clinical situations in order to copy or recreate a prior situation of the patient. This will enable the laboratory technician to overlay the different intraoral scan files of one patient using the CAD software. Insertion, handling and removal of the Di²gital Arch® 1.0 impression screw must be carried out by a qualified professional using the recommended kit and following the guidelines in the catalogues and working protocols.

Check the clinical indications of the Ziacom® Di²gital Arch® 1.0 system thoroughly in the corresponding catalogues and working protocols before planning treatment. Do not use the Ziacom® Di²gital Arch® 1.0 system in locations or for purposes other than the indications described in the catalogues or working protocols.

4. Contraindications

The patient must undergo a preoperative medical examination, as well as assessment of the surgical conditions of the area and of the procedure to determine whether there are any absolute or relative contraindications or risk factors for the use of screws. The user is responsible for evaluating the potential risks and benefits of the use of screws in impression-taking procedures. Ziacom® products should not be used in patients who lack the minimal medical conditions necessary to carry out treatment and prosthetic rehabilitation with dental implants.

General Contraindications:

Patients with systemic diseases (blood clotting disorders, untreated endocrine disorders, rheumatic illnesses, liver disease, neoplasms, diabetes or other serious illness, are taking antiresorptive medicines), psychological or mental illness, poor oral hygiene, severe or uncontrolled periodontal disease, toxic habits (smoking, drug addiction, alcoholism), active infection or pathological processes in the site of placement, unknown allergies to the screw or impression implant materials, nearby anatomical locations (teeth, inferior alveolar canal, nasopalatine duct, mental foramen, nasal sinus, etc.).

Ziacom Implants SLU will update these contraindications as necessary, based on information obtained from clinical studies and post-marketing information.

5. Intended patients

The target population for these medical devices are those patients who, in the opinion of the clinician assessing their case, require treatment in the mandible and maxilla, and are intended to restore chewing function, as well as providing an improvement in occlusion and the distribution of occlusal loads, phonation, aesthetics, prosthetic stability and maintenance of the alveolar process, as well as a positive psychological effect.

These patients must have completed their growth. The age guideline is older than 18 years, this can be confirmed after cessation of growth with a wrist X-ray. Patients at an early age patients will be specifically informed of the changes they may need to make to the prostheses during their lifetime with a with the consequent aesthetic and economic impact due to functional ankylosis of the implant in living tissue.

6. Storage

Ziacom® products must be stored in their original packaging in a dry, clean place away from sunlight and adverse conditions. They must be stored at a temperature below +55°C and above -10°C.

7. Packaging and Supply Format

Ziacom® Di²gital Arch® 1.0 system is supplied in a cardboard box with a product identification label, which functions as a security seal to protect against tampering.

The pack contains:

- The QR code to access these electronic indications for use of the product.

- Plastic heat-sealed blister and labels, indicating the manufacturer, full product reference, batch number, impression screw dimensions and relevant symbols (explained in these instructions for use). The impression screw is inside the blister.

The removable labels are intended for traceability of the screws and/or impression implants, by sticking one label on the patient's history and the remaining one on the clinical record.

Ziacom Implants SLU guarantees that all of its products follow a strict manufacturing, monitoring and cleaning process before being packaged. Before using Ziacom® products, you must check that the packaging is intact and ensure that it is not damaged. If there is any damage or defect, the product must not be used and Ziacom Implants SLU or your authorised official distributor must be notified.

8. Cleaning and Sterilisation

Ziacom® Di²gital Arch® 1.0 impression screws are supplied non-sterile.

Impression screws that are to come into contact with the mouth must be cleaned, disinfected and sterilised before use, following the Ziacom Implants SLU cleaning, disinfection and sterilisation protocol for impression screws, instruments, surgical kit boxes and abutments. You can find this on our website at www.ziacom.com or in Ziacom Implants SLU catalogues.

9. Precautions

Insertion and Instruments

The working protocols describe in detail the precautions and important factors to keep in mind during the impression screw placement process. Prior to placement, it is essential to have consulted the working protocol for the product, along with all other indications regarding the correct use as indicated in the Ziacom® product catalogues.

Special care must be taken to take the necessary steps to prevent the screw, abutments and instruments from being swallowed or aspirated by the patient.

Take extra care not to apply excessive force during placement of the impression screws. If during placement a variation is detected in the bone consistency, stop insertion immediately and make sure that you are not inserting into inadequate bone or dental structure, for example into tooth roots or other anatomical structures, as well as if you detect any deformity in the body of the impression implant that may lead to some type of fracture. It is recommended that a maximum of 50Ncm is not exceeded for Di²gital Arch® 1.0 impression screws.

Attention

It is recommended that the Di²gital Arch® 1.0 impression screws are placed in the mid-vestibular area, both in the maxilla and the mandible. This placement will depend on the clinical condition and the working protocol to be performed.

10. Warnings

Treatment planning for implant and rehabilitation surgery requires specific dental training. The product specifications themselves are no guarantee of proper use. We recommend that users take specialised theoretical and practical training courses in order to understand the proper working techniques and protocols for each product, including the biomechanical, imaging, prosthetic and biological requirements associated with restoration treatment.

It is the user's responsibility to know and find out about the current status of the technique for each product and its possible applications.

Before using Ziacom® products, you must know and be familiar with the corresponding work procedures for placement of Di²gital Arch® 1.0 impression screws, such as their use in scanning procedures with intraoral devices. The patient must meet sufficient anatomical and psychological conditions to undergo treatment with dental implants. The user must ensure proper treatment planning is carried out, to guarantee the safety of the treatment with a minimal margin of error, respecting the patient's vital mouth structures and general health.

The procedure for using the insertion kit for placing Di²gital Arch® 1.0 implants is specified in the corresponding working protocols. The placement of Di²gital Arch® 1.0 impression screws and the digital workflow must be adapted to the individual conditions of the patient, the manufacturer's protocol for the intraoral scanner and the CAD planning software.

Reuse of single-use products carries the risk of possible deterioration of their characteristics. The geometry of the product may have suffered damage from previous use or due to inadequate use. This involves the risk of impression-taking failure and/or damage to the patient's health.

11. Adverse Events or Side Effects

Various adverse or undesired effects can occur during treatment with the Di²gital Arch® 1.0 system. These are documented in the specialised scientific literature published in journals or specialised books on dentistry.

However, the most relevant are:

Lesions in the adjacent anatomical structures (nasal passages, maxillary sinus, vascular structures/nerves), infections, issues with reading the scanbodies and problems inherent to the intraoral scanner and CAD software.

Ziacom Implants SLU will update these side effects as necessary, based on information obtained from clinical studies and post-marketing information.

12. Liability, safety and warranty information

The instructions for use and handling of Ziacom® products are based on published scientific literature, current clinical standards and experience with our products, so they should be considered to be general information intended as a guide. The handling and use of Ziacom® products is the sole responsibility of the user as it is outside the control of Ziacom Implants SLU Ziacom Implants SLU, their subsidiaries and/or their official authorised distributors do not accept any liability, whether explicit or implicit, total or partial, for possible damage or injury caused by improper handling of the product or any other situation not considered in their protocols and instructions for the correct use of their products.

The user must ensure that the Ziacom® product is appropriate for the intended procedure and end purpose. Neither these instructions for use nor the working or handling protocols for the products release the user from this obligation. Ziacom® products must be used, handled and applied clinically by professionals with the appropriate training and qualifications required according to current legislation in each country. Total or partial use, handling or application of any phase of Ziacom®

products by unqualified personnel, or by personnel that does not have the necessary training, will automatically void any type of warranty and may cause serious damage to the patient's health.

Ziacom® products have their own design features and working protocols, so the insertion of the impression screws should be carried out with the adaptors and insertion kit recommended by Ziacom Implants SLU in the catalogue and the working protocols. Insertion of Ziacom® impression screws with adaptors or instruments from other manufacturers may lead to insertion failure, deterioration of the impression screw, adverse impact on the area of insertion in causing damage to the tissues and bone structures. For this reason, insertion should be carried out with the kit and adaptors recommended by Ziacom Implants SLU.

The clinician in charge of the treatment is solely responsible for ensuring the use of original Ziacom® products and that they are used according to the corresponding instructions for use and working protocols throughout treatment. The use of any other non-original Ziacom® instruments, whether alone or in combination with any original Ziacom® products, will immediately void the warranty of the original Ziacom® products. See the Ziacom Implants SLU Warranty Programme (available on the website or by contacting Ziacom Implants SLU, their subsidiaries or authorised distributors).

Any problem or incident with regard to the device should be reported to the manufacturer, Ziacom Implants SLU, to which the user should also send the affected product. In the event of a serious incident, the user should also submit a report to the manufacturer, Ziacom Implants SLU as well as the relevant competent authority, according to the requirements of local legislation.

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.

Warning

Not all Ziacom® products are for sale in every country. Check availability in your country.

The Ziacom® brand and the names of other products and services, including their logos, that are mentioned in this document or on the website www.ziacom.com, are registered trademarks of Ziacom Implants SLU

Ziacom Implants SLU reserves the right to modify, change, eliminate or update any of the products, prices or technical specifications mentioned in this document or in any of its catalogues without prior notification. All rights reserved. The reproduction of this document, whole or in part and in any medium or format, without the corresponding written authorisation from Ziacom Implants SLU is prohibited.

13. Magnetic Resonance

The safety and compatibility of Ziacom® products in a magnetic resonance environment has not been evaluated, nor have they undergone testing to verify whether they heat up or migrate in such environments.

14. Removal of the material

The dental implants, abutments and kit, as well as each product and consumable used during surgery for placing dental implants and which are disposed of may pose a health hazard to people handling them. For this reason, you are advised to check and comply with current legislation on biohazardous waste.

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EN | CLEANING, DISINFECTION AND STERILISATION

Cleaning, disinfection and sterilisation

The protocols described in this section must only be carried out by personnel qualified to clean, disinfect and sterilise the dental materials specified here in.

CLEANING AND DISINFECTION INSTRUCTIONS

Applicable for instruments, surgical and prosthetic boxes.

Disassembly

1. Dismount* the appropriate instruments, for example manual ratchets, drills or drill stops.
2. Remove the various components from the surgical or prosthetic box for correct cleaning.

Cleaning and disinfection

For disinfecting instruments and surgical boxes:

1. Submerge the instruments in a detergent/disinfectant solution** suitable for dental instruments to help eliminate any adhered biological residues. If an ultrasound bath is available***, confirm that the detergent/disinfectant solution is indicated for use with this type of equipment.
2. Manually remove any biological residues with a non-metallic brush and pH-neutral detergent.
3. Rinse with copious water.
4. When cleaning the surgical and prosthetic boxes, always use a pH-neutral detergent and non-abrasive utensils to avoid damaging the surface of the boxes.
5. Dry the materials with disposable cellulose, lint-free clothes or compressed air.

For disinfecting plastic caps and spacers:

1. Submerge in a neat benzalkonium chloride solution for 10 minutes.
2. Rinse with distilled water.
3. Dry the caps and spacer before use.

Inspection

1. Check that the instruments are perfectly clean; if not, repeat the cleaning and disinfection steps.
2. Discard any instruments with imperfections and replace them before the next procedure.
3. Check that the instruments and the surgical and prosthetic boxes are perfectly dry before reassembling the parts and proceeding to their sterilisation.

* See assembly and disassembly instructions.

** Follow the instructions from the disinfectant's manufacturer to determine the correct concentrations and times.

*** Follow the instructions from the ultrasound bath's manufacturer to determine the correct temperature, concentration and times.

STERILISATION INSTRUCTIONS FOR STEAM AUTOCLAVE

Applicable to orthodontic implants, abutments, instruments, surgical and prosthetic boxes, pins, fixation screws and laminates.

1. Introduce each material separately in individual sterilisation bags, then seal the bags. For joint sterilisation, place the instruments in their surgical box, introduce the box into a sterilisation bag and seal the bag.
2. Place the bags to be sterilised in the autoclave.
3. Sterilise in a steam autoclave at 134°C/273°F (max. 137°C/276°F) for 4 min (minimum) and at 2 atm. Torque wrenches must be sterilised in 3 vacuum cycles at 132°C/270°F for a minimum of 1.5 minutes and vacuum-dried for a minimum of 20 minutes.

For the United States only: The validated and recommended sterilisation cycle for the US must be performed in a steam autoclave at 132°C/270°F for at least 15 min and with the drying time of at least 15- 30 min.

Important:

- Make sure the drying stage is allowed to run to completion, otherwise the products may be damp.
- Check the sterilisation equipment if the materials or sterilisation bags are damp at the end of the sterilisation cycle.
- Perform the necessary maintenance actions on the autoclave according to the established periodicity and following the manufacturer's instructions.

STORAGE OF ZIACOM® PRODUCTS

- Store the products in their original packaging and in a clean, dry location until they are used.
- After sterilisation, keep the products in the sealed sterilisation bags and in a clean, dry location.
- Never exceed the use by date indicated by the manufacturer of the sterilisation bags.
- Always follow the indications of the manufacturer of the sterilisation bags.

GENERAL RECOMMENDATIONS

- Never use damaged or dirty material; never reuse single-use products. The user is responsible for following the instructions described in this document correctly.
- The attention to piercing or sharp elements. Gloves should be worn when cleaning the materials to avoid accidents during handling.
- Follow the safety instructions indicated by the manufacturer of the disinfectant agent.
- The product's sterility cannot be guaranteed if the sterilisation bag is open, damaged or damp.
- Respect all stages of the sterilisation process. If the materials or sterilisation bags contain traces of water or moisture, check the autoclave and repeat the sterilisation.
- Orthodontic abutments and implants are supplied UNSTERILISED and must always be sterilised before use.
- Instruments and surgical and prosthetic boxes are supplied UNSTERILISED and must always be sterilised before use and cleaned and disinfected after use.
- The sterilisation, cleaning and disinfection processes gradually deteriorate the instruments. Inspect the instruments thoroughly to detect any signs of deterioration.
- Avoid contact between products made from different materials (steel, titanium, etc.) during the cleaning, disinfection and sterilisation processes.
- Ziacom Implants SLU recommends these instructions are implemented for the correct maintenance and safety of their products; accordingly, the company refuses any liability for any damage to the products that could arise if the user applies alternative cleaning, disinfection and sterilisation procedures.

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PRESENTATION AND HANDLING OF THE DI2GITAL ARCH® 1.0 FIXING SCREW.

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EN Press the word “PRESS” and open the carton. Remove the flap from the carton and take out the blister pack.

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EN Carefully remove the seal from the blister pack.

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
EN Drop the Di2gital Arch® 1.0 impression screw onto a clean cloth.

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EN Sterilise the Di2gital Arch® 1.0 screws and the kit according to the cleaning, disinfection and sterilisation instructions.

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
EN Once the screws and instruments are sterile, use the manual or contra-angle insertion tool to hold the Di2gital Arch® 1.0 fixing screw for clinical use.

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EN Remember to remove the identification labels from the Di2gital Arch® 1.0 impression screw to stick them onto the patient's clinical records, thereby enabling traceability of the product.

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Ziacom Medical USA LLC