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EN | INSTRUCTIONS FOR USE

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FR | INSTRUCTIONS D'USAGE

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IMPLANT INSTRUCTIONS
FOR USE

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Ziacom[®] Di²gital Arch[®]

INSTRUCTIONS FOR USE

Instructions for use applicable to the implant systems Di²gital Arch[®] 2.0.

IMPORTANT INFORMATION

Please read these instructions carefully before using the product.

This document contains basic information on the use of the Ziacom[®] original Di²gital Arch[®] 2.0 implant, hereinafter referred to as the Ziacom Di²gital Arch[®] system, or simply "Ziacom[®] products". This documentation has been prepared as a quick guide to be consulted by the clinician responsible for treatment, hereinafter the "user", and does not provide sufficient indications or technical specifications for correct use of Ziacom[®] products. This is not an alternative nor substitute for specialised training or professional clinical experience.

Ziacom[®] products must be used according to a suitable treatment plan and in strict compliance with the manufacturer's working protocols: Ziacom Medical S.L. Please read the specific working protocols for this product, as well as the instructions for use and cleaning, disinfection and sterilisation carefully before using this Ziacom[®] product. You can find this information on our website at www.ziacom.com or request it from your Ziacom Medical S.L. authorised distributor.

1. Description of the System

Ziacom[®] products consist of a Di²gital Arch[®] the 2.0 implant, digital impression abutments and surgical kit.

Ziacom Medical S.L. has created a Di²gital Arch[®] 2.0 implant, which serves as a reference structure to facilitate alignment of virtual files in various clinical situations in which you wish to copy the prior situation of the patient in prosthetic procedures with intraoral scanners, in accordance with the existing scientific literature and current clinical standards. Ziacom Medical S.L. also develops, manufactures and markets latest-generation implant surgery products and a wide range of restoration solutions to meet each patient's particular needs.

2. Description of the Products

Di²gital Arch[®] 2.0 System

The Ziacom[®] Di²gital Arch[®] 2.0 implant is manufactured in biocompatible grade-4 commercially pure titanium, in accordance with standards ASTM F67-13 and UNE EN-ISO 5832-2:2018, with surface treatment that speeds up osseointegration of the implant.

3. Indications for Use

Before carrying out treatment with Ziacom[®] products, it is vital to assess the patient beforehand with clinical diagnosis, imaging and study models, as well as determine possible contraindications for the treatment. Correct impression taking involves

creating a prosthesis and planning the restoration treatment, addressing functional defects (phonation, mastication etc.).

Di2gital Arch® 2.0 System

The Ziacom® Di²gital Arch® 2.0 system is indicated to serve as a reference structure when aligning the virtual files of a single patient at different stages of treatment or in different clinical situations in order to copy or recreate a prior situation of the patient. This will enable the laboratory technician to overlay the different intraoral scan files of one patient using the CAD software. Insertion and handling of the Di²gital Arch® 2.0 implant must be carried out by a qualified professional using the recommended kit and following the guidelines in the catalogues and working protocols.

Important

Check the clinical indications of the Ziacom® Di²gital Arch® 2.0 system thoroughly in the corresponding catalogues and working protocols before planning treatment. Do not use the Ziacom® Di²gital Arch® 2.0 system in locations or for purposes other than as described in the indications in the catalogues or working protocols.

4. Contraindications

The patient must undergo a preoperative medical examination, as well as assessment of the surgical conditions of the area and of the procedure to determine whether there are any absolute or relative contraindications or risk factors for the use of implants. The user is responsible for evaluating the potential risks and benefits of the use of the implants in impression-taking procedures. Ziacom® products should not be used in patients who lack the minimal medical conditions necessary to carry out treatment and prosthetic rehabilitation with dental implants.

General Contraindications:

Patients with systemic diseases (blood clotting disorders, untreated endocrine disorders, rheumatic illnesses, liver disease, neoplasms, diabetes or other serious illness, are taking antiresorptive medicines), psychological or mental illness, poor oral hygiene, severe or uncontrolled periodontal disease, toxic habits (smoking, drug addiction, alcoholism), active infection or pathological processes in the site of placement, unknown allergies to the screw or impression implant materials, nearby anatomical locations (teeth, inferior alveolar canal, nasopalatine duct, mental foramen, nasal sinus, etc.).

Ziacom Medical S.L. will update these contraindications as necessary, based on information obtained from clinical studies and post-marketing information.

5. Storage

Ziacom® products must be stored in their original packaging in a dry, clean place away from sunlight and adverse conditions. They must be stored at a temperature below +55°C and above -10°C.

6. Packaging and Supply Format

Ziacom® Di²gital Arch® 2.0 system is supplied in a cardboard box with a product identification label, which functions as a security seal to protect against tampering.

The pack contains:

- The QR code to access these electronic indications for use of the product.
- Plastic heat-sealed blister and labels, indicating the manufacturer, full product reference, batch number, impression screw dimensions and relevant symbols (explained in these instructions for use). The implant is inside the blister.

The removable labels are intended for traceability of the impression implants, by sticking one label on the patient's history and the remaining one on the clinical record.

Ziacom Medical S.L. guarantees that all of its products follow a strict manufacturing, monitoring and cleaning process before being packaged. Before using Ziacom® products, you must check that the packaging is intact and ensure that it is not damaged. If there is any damage or defect, the product must not be used and Ziacom Medical S.L. or your authorised official distributor must be notified.

7. Cleaning and Sterilisation

The Ziacom® Di²gital Arch® 2.0 is supplied sterilised using 25 kGy beta-radiation sterilisation. If the sterile container is altered, tampered with, damaged or broken, sterility is not guaranteed and the product must not be used. Do not use the product after the expiry date shown on the container.

Resterilisation of Ziacom® dental implants is NOT authorised.

Important

The outside of blister pack containing the implant vial is NOT sterile; only the inside of the blister pack is sterile. Do not place the blister pack in the sterile surgical field. It must NOT be opened until just before placement of the implant.

8. Precautions

Insertion and Instruments

The working protocols describe in detail the precautions and important factors to keep in mind during the process for placing impression implants. Prior to placement, it is essential to have consulted the working protocol for the product, along with all other indications regarding the correct use as indicated in the Ziacom® product catalogues.

Special care must be taken to take the necessary steps to prevent the implant, abutments and instruments from being swallowed or aspirated by the patient.

Drilling must always be performed by controlling the axis of drilling, tapping and under irrigation and at the revolutions per minute indicated in the working protocol.

Take extra care not to apply excessive force during placement of the impression implants. If during placement a variation is detected in the bone consistency, stop insertion immediately and make sure that you are not inserting into inadequate bone

or dental structure, for example into tooth roots or other anatomical structures, as well as if you detect any deformity in the body of the impression implant that may lead to some type of fracture. It is recommended that a maximum of 50Ncm is not exceeded for Di²gital Arch® 2.0 implants.

Attention

It is recommended that the Di²gital Arch® 2.0 implants are placed in the mid-vestibular area, both in the maxilla and the mandible. This placement will depend on the clinical condition and the working protocol to be performed.

9. Warnings

Treatment planning for implant and rehabilitation surgery requires specific dental training. The product specifications themselves are no guarantee of proper use. We recommend that users take specialised theoretical and practical training courses in order to understand the proper working techniques and protocols for each product, including the biomechanical, imaging, prosthetic and biological requirements associated with restoration treatment.

It is the user's responsibility to know and find out about the current status of the technique for each product and its possible applications.

Before using Ziacom® products, you must know and be familiar with the corresponding work procedures for placement of Di²gital Arch® 2.0 implants, such as their use in scanning procedures with intraoral devices. The patient must meet sufficient anatomical and psychological conditions to undergo treatment with dental implants. The user must ensure proper treatment planning is carried out, to guarantee the safety of the treatment with a minimal margin of error, respecting the patient's vital mouth structures and general health.

The procedure for using the insertion kit for placing Di²gital Arch® 2.0 implants is specified in the corresponding working protocols. The placement of Di²gital Arch® 2.0 implants and the digital workflow must be adapted to the individual conditions of the patient, the manufacturer's protocol for the intraoral scanner and the CAD planning software.

Reuse of single-use products carries the risk of possible deterioration of their characteristics. The geometry of the product may have suffered damage from previous use or due to inadequate use. This involves the risk of impression-taking failure and/or damage to the patient's health.

10. Adverse Events or Side Effects

Various adverse or undesired effects can occur during treatment with the Di²gital Arch® 2.0 system. These are documented in the specialised scientific literature published in journals or specialised books on dentistry.

However, the most relevant are:

Lesions in the adjacent anatomical structures (nasal passages, maxillary sinus, vascular structures/nerves), infections, failure of osseointegration of the Di²gital Arch® 2.0 implant, issues with reading the scanbodies and problems inherent to the intraoral scanner and CAD software.

Ziacom Medical S.L. will update these side effects as necessary, based on information obtained from clinical studies and post-marketing information.

11. Liability, safety and warranty information

The instructions for use and handling of Ziacom® products are based on published scientific literature, current clinical standards and experience with our products, so they should be considered to be general information intended as a guide. The handling and use of Ziacom® products is the sole responsibility of the user as it is outside the control of Ziacom Medical S.L. Ziacom Medical S.L., their subsidiaries and/or their official authorised distributors do not accept any liability, whether explicit or implicit, total or partial, for possible damage or injury caused by improper handling of the product or any other situation not considered in their protocols and instructions for the correct use of their products.

The user must ensure that the Ziacom® product is appropriate for the intended procedure and end purpose. Neither these instructions for use nor the working or handling protocols for the products release the user from this obligation. Ziacom® products must be used, handled and applied clinically by professionals with the appropriate training and qualifications required according to current legislation in each country. Total or partial use, handling or application of any phase of Ziacom® products by unqualified personnel, or by personnel that does not have the necessary training, will automatically void any type of warranty and may cause serious damage to the patient's health.

Ziacom® products have their own design features and working protocols, so the insertion of the impression implants should be carried out with the adaptors and insertion kit recommended by Ziacom Medical S.L. in the catalogue and the working protocols. Insertion of Ziacom® impression implants with adaptors or instruments from other manufacturers may lead to insertion failure, deterioration of the impression implant, adverse impact on the area of insertion in causing damage to the tissues and bone structures. For this reason, insertion should be carried out with the kit and adaptors recommended by Ziacom Medical, S.L.

The clinician in charge of the treatment is solely responsible for ensuring the use of original Ziacom® products and that they are used according to the corresponding instructions for use and working protocols throughout treatment. The use of any other non-original Ziacom® instruments, whether alone or in combination with any original Ziacom® products, will immediately void the warranty of the original Ziacom® products. See the Ziacom Medical S.L. Warranty Programme (available on the website or by contacting Ziacom Medical SL, their subsidiaries or authorised distributors).

Warning

Not all Ziacom® products are for sale in every country. Check availability in your country.

Ziacom®, Basic®, Di2gital Arch®, DSQ®, Galaxy®, Kiran®, Kirator®, OsseosBCP®, OsseosTCP®, PlexGuide®, Titansure®, Titansure® Active, Tx30®, XDrive®, Zellplex®, ZiaCam®, Ziacom®3D, Ziacor®, Zinic®, Zinic®MT, Zinic®Shorty, ZiTium®, ZM1®, ZM4®, ZM4®MT, ZM8®, ZM8®N, ZM8®S, ZM-Equator®, ZMK®, ZMR®, ZPlus®, ZV2® and all their logos are registered trademarks of Ziacom Medical S.L.

The Ziacom® brand is a registered trademark, as are some of the products that may or may not be mentioned in this document. Ziacom Medical S.L. reserves the right to modify, change, eliminate or update any of the products, prices or technical specifications mentioned in this document or in any of its catalogues without prior notification. All rights reserved. The reproduction of this document, whole or in part and in any medium or format, without the corresponding written authorisation from Ziacom Medical SL is prohibited.

12. Magnetic Resonance

The safety and compatibility of Ziacom® products in a magnetic resonance environment has not been evaluated, nor have they undergone testing to verify whether they heat up or migrate in such environments.

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Ziacom®, Basic®, Di2igital Arch®, DSQ®, Galaxy®, Kiran®, Kirator®, OsseosBCP®, OsseosTCP®, PlexGuide®, Titansure®, Titansure® Active, Tx30®, XDrive®, Zellplex®, ZiaCam®, Ziacom®3D, Ziacor®, Zinic®, Zinic®MT, Zinic®Shorty, ZiTium®, ZM1®, ZM4®, ZM4®MT, ZM8®, ZM8®N, ZM8®S, ZM-Equator®, ZMK®, ZMR®, ZPlus®, ZV2® e todos os seus logótipos são marcas registadas da Ziacom Medical SL.

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Symbols of labelling codes and used in catalogues, standard ISO 15223-1

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Rx Only

Caution: federal law prohibite dispensing without prescription.

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


EN Take out the microimplant set – Di2gital Arch® 2.0 insertion tool, pulling upwards vertically on the vial for its clinical use.

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EN Remember to remove the identification labels from the Di2gital Arch® 2.0 fixing screw to stick them onto the patient's clinical records, thereby enabling traceability of the product.

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